Simplifying Proteomics

PTMScan® Products and Services

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
- 3 Key to Success: PTMScan® Antibodies
- 4 Discovery vs. Direct
- 5 PTMScan® Kits vs. Services
- 6 PTMScan® Discovery and Direct Services
- 7 Complementary Services
- 8 PTMScan® Discovery Products and Services
- 13 PTMScan® Direct Services
- 14 From Bench to Bedside
- 15 Support for Drug Discovery

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



Serum-starved H₂O₂ Treated H₂O₂ + λ phosphatase

Phospho-Tyrosine (P-Tyr-1000) MultiMab™ Rabbit mAb mix #8954: Confocal IF analysis of C2C12 cells, serum-starved (left), treated with H₂O₂ (2 mM, 10 min; middle), or treated with H₂O₂ followed by λ phosphatase (right), using #8954 (green). Actin filaments were labeled with DyLight™ 554 Phalloidin #13054 (red). Blue pseudocolor = DRAQ5® #4084 (fluorescent DNA dye).

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
PTMScan® Discovery Products and Services

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UbiScan® – Ubiquitination Proteomics

Target	Motif	Antibody [‡]	PTMScan®		KinomeView®	
		No.	Kit	Service	Kit #9812	Service
Ubiquitin Remnant	K-ε-GG	–	#14482/#5562	✓	–	✓

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Sook-Young, J. et al. (2015) Proteomic Analysis of Ubiquitin-Like Post-translational Modifications Induced by the Adenovirus E4-ORF3 Protein. *Journal of Virology* 89(3) 1744-55.

Kronke, J. et al. (2015) Lenalidomide induces ubiquitination and degradation of CK1α in del(5q) MDS. *Nature* 523(75599), 183–188.

Kronke, J. et al. (2014) Lenalidomide causes selective degradation of IKZF1 and IKZF3 in multiple myeloma cells. *Science* 343(6168), 301–305.

Theurillat, J.P. et al. (2014) Ubiquitylome analysis identifies dysregulation of effector substrates in SPOP-mutant prostate cancer. *Science* 346(6205), 85–89.

Mathew, R. et al. (2012) BTB-ZF factors recruit the E3 ligase cullin 3 to regulate lymphoid effector programs. *Nature* 491(7425), 618–621.

Emanuele, M.J. et al. (2011) Global identification of modular Cullin-RING ligase substrates. *Cell* 147(2), 459–474.

Kim, W. et al. (2011) Systematic and quantitative assessment of the ubiquitin-modified proteome. *Mol. Cell* 44(2), 325–340.

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

Area proportional Venn diagram showing the overlap between motif antibody and IMAC enrichment of phosphopeptides from MKN-45 cells (A). Comparisons were made on unique proteins/sites. MKN-45 cells were treated and profiled as outlined (B). Combinations of motif antibodies were used to enrich phosphopeptides as indicated (C).

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Paardekooper Overman, J. et al. (2014) PZR Coordinates Shp2 Noonan and LEOPARD syndrome signaling in zebrafish and mice. *Mol. Cell Biol.* 34(15), 2874–2889.

Siddoway, B. et al. (2013) Synaptic activity bidirectionally regulates a novel sequence-specific S-Q phosphoproteome in neurons. *Journal of Neurochemistry* 128(6):841–851.

Ren, H., et al. (2012) Identification of anaplastic lymphoma kinase as a potential therapeutic target in ovarian cancer. *Cancer Res.* 72(13) 3312–3323.

Carretero, J. et al. (2010) Integrative genomic and proteomic analyses identify targets for Lkb1-deficient metastatic lung tumors. *Cancer Cell* 17(6), 547–559.

Moritz, A. et al. (2010) Akt-RSK-S6 kinase signaling networks activated by oncogenic receptor tyrosine kinases. *Sci. Signal* 3(136), ra64.

Rikova, K. et al. (2007) Global survey of phosphotyrosine signaling identifies oncogenic kinases in lung cancer. *Cell* 131(6), 1190–1203.

Griswold, I.J., et al. (2006) Kinase domain mutants of Bcr-Abl exhibit altered transformation potency, kinase activity, and substrate utilization, irrespective of sensitivity to imatinib. *Mol Cell Biol.* 26(16) 6082–6093.

Rush, J. et al. (2005) Immunoaffinity profiling of tyrosine phosphorylation in cancer cells. *Nat. Biotechnol.* 23(1), 94–101.

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


Western blot analysis of lysates from untreated (-) or AdOx-treated (+) HCT116 cells using asymmetric dimethyl arginine (R-2Me(a)) antibodies. The specific signal was blocked by antigen R-2Me(a) peptides (middle panel), but not by corresponding symmetric dimethyl arginine (R-2Me(s)) peptides (right panel).*

*This research was originally published in Molecular and Cellular Proteomics. Guo, A., et al. Immunoaffinity enrichment and mass spectrometry analysis of protein methylation. *Mol. Cell. Biol.* 2014; 13(1):372-387. © the American Society for Biochemistry and Molecular Biology, and is licensed under CC BY 4.0 <http://creativecommons.org/licenses/by/4.0/>.

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MethylScan® – Methylation Proteomics

Target	Motif	Antibody [‡]	PTMScan®		KinomeView®	
		No.	Kit	Service	Kit #9812	Service
Mono-Methyl Arginine	R-Me	#8015	#12235	✓	—	✓
Asymmetric Di-Methyl Arginine	R-2Me(a)	#13522	#13474	✓	—	✓
Symmetric Di-Methyl Arginine	R-2Me(s)	#13222	#13563	✓	—	✓
Pan-Methyl Lysine	K-Me/K-2Me/K-3Me	#14117 (2Me) #14680 (3Me)	#14809	✓	—	✓

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MethylScan®の使用文献

Guo, A. et al. (2014) Immunoaffinity enrichment and mass spectrometry analysis of protein methylation. *Mol. Cell. Proteomics* 13(1) 372-387.

Sylvestersen, K. B. et al. (2014) Proteomic analysis of arginine methylation sites in human cells reveals dynamic regulation during transcriptional arrest. *Mol. Cell. Proteomics* 13(8):2072-2088.

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

Target	Motif	Antibody [‡]		PTMScan®		KinomeView®	
		No.	Kit	Service	Kit #9812	Service	
AcetylScan® – Acetylation Proteomics							
Acetyl-Lysine	Ac-K	#9814	#14499/#13416	✓	–	–	✓
GlutarylScan™ – Glutarylation Proteomics							
Glutaryl-Lysine	Glut-K	–	#26101	✓	–	–	✓
MalonylScan™ – Malonylation Proteomics							
Malonyl-Lysine	Mal-K	–	#93872	✓	–	–	✓
PropionylScan™ – Propionylation Proteomics							
Propionyl-Lysine	Prop-K	–	#17848	✓	–	–	✓
SuccinylScan™ – Succinylation Proteomics							
Succinyl-Lysine	Succ-K	–	#13764	✓	–	–	✓

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Svinkina, T. et al. (2015) Deep, quantitative coverage of the lysine acetylome using novel anti-acetyl-lysine antibodies and an optimized proteomic workflow. *Mol. Cell Proteomics* 14(9), 2429–2440.

Bouchut, A. et al. (2015) Proteome-wide lysine acetylation in cortical astrocytes and alterations that occur during infection with brain parasite *Toxoplasma gondii*. *PLoS ONE* 10(3), e0117966.

Ryder, D.J. et al. (2015) Identification of the acetylation and ubiquitin-modified proteome during the progression of skeletal muscle atrophy. *PLoS ONE* 10(8), e0136247.

Rardin, M.J. et al. (2013) SIRT5 regulates the mitochondrial lysine succinylome and metabolic networks. *Cell Metab.* 18(6), 920–933.

Mielcarek, M. et al. (2013) HDAC4 Does Not Act as a Protein Deacetylase in the Postnatal Murine Brain In Vivo. *PLoS One* 8(11), e80849.

Jeffers, V. et al. (2012) Lysine acetylation is widespread on proteins of diverse function and localization in the protozoan parasite *Toxoplasma gondii*. *Eukaryot. Cell* 11(6) 735–742.

PTMScan® Discovery Products and Services (cont.)

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Caspase Cleavage Substrate Proteomics		Antibody [‡]	PTMScan®		KinomeView®	
Target	Motif	No.	Kit	Service	Kit #9812	Service
Caspase Cleavage Substrate	[DE(T/S/A/D)]	#8698	#12810	✓	-	✓

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Anania, V.G. and Lill, J.R. (2015) Proteomic tools for the characterization of cell death mechanisms in drug discovery. *PROTEOMICS - Clinical Applications* 9(7-8) 671-683.

Pham, V.C. et al. (2012) Complementary proteomic tools for the dissection of apoptotic proteolysis events. *J. Proteome Res.* 11(5), 2947-2954.

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

- Adherens Junction Dynamics
- AMPK Signaling
- B Cell Receptor Signaling
- Cell Cycle Control: G1/S Checkpoint
- Cell Cycle Control: G2/M DNA Damage Checkpoint
- ErbB/HER Signaling
- G Protein-coupled Receptor Signaling to MAPK/Erk
- Insulin Receptor Signaling
- Jak/Stat Signaling: IL-6 Receptor Family
- MAPK Signaling Cascades
- MAPK/Erk in Growth and Differentiation
- NF-κB Signaling
- PI3 Kinase/Akt Signaling*
- Regulation of Actin Dynamics
- Regulation of Microtubule Dynamics
- SAPK/JNK Signaling Cascades
- Signaling Pathways Activating p38 MAPK
- T cell Receptor Signaling
- Toll-like Receptor Signaling

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PTMScan® Directの使用文献

Stokes, M.P. et al. (2013) Quantitative profiling of DNA damage and apoptotic pathways in UV damaged cells using PTMScan Direct. *Int. J. Mol. Sci.* 14(1), 286–307.

Pease, B.N. et al. (2013) Global Analysis of Protein Expression and Phosphorylation of Three Stages of Plasmodium falciparum Intraerythrocytic Development. *J. Proteome Res.* 2(9):4028–4045.

Stokes, M.P. et al. (2012) PTMScan Direct: identification and quantification of peptides from critical signaling proteins by immunoaffinity enrichment coupled with LC-MS/MS. *Mol. Cell Proteomics* 11(5), 187–201.

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


ALK (D5F3®) XP® Rabbit mAb #3633: IHC analysis of paraffin-embedded human lung carcinoma with high (left) and low levels (right) of ALK fusion protein expression using #3633.

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1. Rikova, K. et al (2007) *Cell* 131(6), 1190-203.
2. Soda, M., et al (2007) *Nature* 448(7153), 561-6.
3. Koivunen, J.P., et al (2008) *Clin. Cancer Res.* 14(13), 4275-83.
4. Shaw, A.T., et al (2009) *J. Clin. Oncol.* 27(26), 4247-53.
5. Press Release: <http://www.fda.gov/NewsEvents/Newsroom/PressAnnouncements/ucm269856.htm>
6. Press Release: <http://www.ventana.com/site/page?view=press-release-jan21-2015>

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Troy, Product Scientist,
has been with CST since 2010.

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www.cstj.co.jp/support

© 2015 Cell Signaling Technology, Inc.
Cell Signaling Technology, CST, AcetylScan,
GlutarylScan, KinomeView, MalonoylScan,
MethylScan, MultiMab, PathScan, PhosphoScan,
PropionylScan PTMScan, SignalSilence,
SuccinylScan, UbiScan, XP, and D5F3 are
trademarks of Cell Signaling Technology, Inc.
All other trademarks are the property of their
respective owners.

♻️ Printed on FSC® certified paper using vegetable inks and processed chlorine free.

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www.cstj.co.jp

CST Antibody Performance Guarantee: To learn more, please visit: www.cellsignal.com/abguarantee.

COVER IMAGE: Cellular Landscape: Vesicle Trafficking www.cellsignal.com/VesicleTrafficking.

16BROPTMS0136JPN For Research Use Only. Not For Use in Diagnostic Procedures.